

Pilot Core Drilling Machine

RDC-179

Electrodes

Lining

R&D

General Description

The evaluation of raw materials, green pastes and electrodes under laboratory conditions is advantageous in terms of time consumption and cost optimizations. The pilot electrodes can be prepared with the RDC-160 and RDC-165 equipment. To obtain samples for the measurement of their properties, the baked pilot electrodes need to be core drilled, the aim being to obtain samples with a 50 mm diameter and a length of approximately 200 mm.

The core drilling is conducted with the RDC-179 apparatus, where the pilot electrode is first clamped. A core bit, driven by an electrical motor, automatically drills a core sample with a given speed, while water is sprayed to ensure a smooth core geometry and surface. No further operation from the operator is required, except for placing the electrode in the right position for drilling.

After drilling, the cores are ready for the cutting operation (see equipment RDC-140, RDC-148 and RDC-149), depending on the properties to analyze.



Technical Information

Function	Drilling of Cores in Pilot Electrodes
Sample	Pilot electrode up to Ø 180 mm x 200 mm
Process Time	~ 5 minutes
Installation	Floor standing
Dimensions (LxWxH)	140 x 80 x 316 cm
Weight	465 kg
Electrical Property	400 V 3/N/PE, 50 Hz 2.2 kW, 5.5 A
Fluid Property	Water : 480 l/h, 3-7 bar
Database Connection	No

Additional Recommended Equipment

Saw (RDC-140 or RDC-148 or RDC-149)
Drying oven (min. temperature 180°C)
Test equipment for analysis