

Bench Drilling Machine

RDC-180

Electrodes

Lining

Lab

General Description

The anode air reactivity is measured with the RDC-151 equipment. For this test, a center hole of a given depth is required in the face side of the sample.

The center hole can be drilled with the RDC-180 apparatus. It consists of a drilling machine specifically adjusted to meet the geometry requirements of the air reactivity test.



Technical Information

Function	Sample Preparation for RDC-151
Sample	Core Ø 50 x 60 mm
Process Time	~ 1 minutes
Installation	Workbench
Dimensions (LxWxH)	23 x 39 x 84 cm
Weight	27 kg
Electrical Property	230 V 1/N/PE, 50 Hz 0.54 kW, 2.5 A
Database Connection	No

Additional Recommended Equipment

Air Reactivity Anode (RDC-151)
Vacuum cleaner